

This material is provided for educational use only. The information in these slides including all data, images and related materials are the property of :

Andrew Ward and Liz Kubalek

Department of Cell Biology

The Scripps Research Institute

10550 North Torrey Pines Rd

La Jolla, CA 92037, USA



Tel: (858) 784- 9827

Fax: (858) 784- 2749

No part of this material may be reproduced without explicit written permission.


Crystallization of Membrane Proteins

Reconstitution






GalCer Scaffolding

- Create tubes made of bilayers doped with a tag (Ni^{2+} , charge, biotin, etc.)
- Add tagged protein and lipid in the presence of detergent to pre-formed tubes
- Dialyze detergent out to form a second bilayer made of a reconstituted helical array of protein



GalCer Experiments



no dialysis

with dialysis

Detergent Removal Strategies




- Dialysis membranes
 - slow, difficult to set up with small volumes
- Dilution
 - very fast, can change [x] unfavorably
- BioBeads
 - fast, easy, can vary rate many ways


BioBead preparation

- Clean BioBeads with methanol
- Wash BioBeads in buffer (+/- detergent)
- Store BioBeads wet on moist filter paper for up to one month


Detergent Removal via BB



2 hours



6 hours




16 hours

Screens




Set up basic matrices to begin

- pH, [lip]:[prot], lipid choice
- Refine with additives, variations in [x], detergent removal rate, temp, etc.


Crystallogenesis



Crystallogenesis




Helical Arrays



GalCer tubes:


- confined geometry
- susceptible to detergent
- template for any type helix



Native tubes:

- protein must drive helices
- slower?


Helical Crystals




2D Arrays



2D Crystals



2D Crystals



What makes a good crystal?

- Low mosaicity/no variation in helical pitch
- Larger (2D)/longer, straighter (helical)
- Low background
- System amenable to cryo...especially for helical arrays